

The new MedSurg Pro® by DARCO

New design – high level of functionality



- > postoperative for stabilisation and off-loading
- > right and left fit
- > high and low version available

Toe cap optional available



DARCO

MedSurg Pro[®]



- > extending the ankle
- > optimized heel guidance
- > large heel counter for improved rear foot stability

> light, breathable u

The MedSurg Pro[®] is the optimized successor model of the MedSurg healing shoe by DARCO. The shoe provides stabilisation and off-loading after surgery – now available in two versions.

The light and breathable upper together with the smooth and soft lining allow maximum comfort. Strap loops enable an easy use of the velcro[®] closure system providing high stability.

Featuring quality classics as square-toe design, low profile, wide entry and distal joint spring (metatarsal-metashank), the shoe maintains its proven high standard on a substantial basis.



> non-slippery sole

> removable insole
> exchangeable with
or Puzzle Insole[®]





upper



- > broad velcro® fastener across the instep
- > strap loops on both sides of the shoe



- > wide entry
- > easy use



with PegAssist® Insole

- > distal joint spring for rigidity in the forefoot area
- > more than 10 cm sole length stability

- > extra space for removable inserts (also diabetic proof foot beddings)
- > low profile – no height adjustment necessary



Article List

Item-No.	Colour	Size
----------	--------	------

High version:

Women:

MS-W1H	Black	WS	34.0 – 36.5
MS-W2H	Black	WM	37.0 – 38.5
MS-W3H	Black	WL	39.0 – 41.0

Men:

MS-M1H	Black	MS	39.0 – 41.0
MS-M2H	Black	MM	41.5 – 43.0
MS-M3H	Black	ML	43.5 – 45.0
MS-M4H	Black	MXL	45.5 – 47.0

Low version:

Women:

MS-W1	Black	WS	34.0 – 36.5
MS-W2	Black	WM	37.0 – 38.5
MS-W3	Black	WL	39.0 – 41.0

Men:

MS-M1	Black	MS	39.0 – 41.0
MS-M2	Black	MM	41.5 – 43.0
MS-M3	Black	ML	43.5 – 45.0
MS-M4	Black	MXL	45.5 – 47.0

Sizes and Measurements



WS 34–36.5




WM 37–38.5



WL 39–41



MS 39–41



MM 41.5–43



ML 43.5–45



MXL 45.5–47

All dimensions are shown in cm and relate to the insoles.

